

EYE TEST CHART

Remove your glasses and step back approximately 14 inches away from the eye test chart. Read the lines starting at the top, the first line that you have trouble reading has the lens' strength to the left. That is your correct strength. All glasses on this display have the strength indicated on them and are color-coded to make finding your strength easier.

+3.25	The Clear Choice
+2.75	The Clear Choice
+2.50	The Clear Choice
+2.25	The Clear Choice
+2.00	The Clear Choice
+1.75	The Clear Choice
+1.50	The Clear Choice
+1.25	The Clear Choice
+1.00	The Clear Choice

Attention: Ready-to-wear non prescriptive glasses are not intended to replace prescribed corrective lenses or examinations by an eye care professional. Continuous eye check-ups are necessary to determine your eye health status and vision needs.

Print eye chart at actual size